

# CleanShield® Endocavity Ultrasound Storage Cabinet

The CleanShield® Ultrasound Probe Storage Cabinet is designed to effectively and securely store disinfected endocavity transducers and semi-critical devices by suspending each disinfected probe in a vertical position that minimizes stress on the connection cable, strain relief and electrical pack of the probe. Specifically designed for vaginal, rectal and general purpose ultrasound probes.



Cat. # CS-VR06

## SPECIFICATIONS

	CS-VR06
Width	26"
Depth	10"
Height	88"

## SPECIFICATIONS

	CS-VR06	CS-VR06 220
<b>Construction</b>		
Main Chamber	Polypropylene	Polypropylene
Door	Polypropylene with polycarbonate window	Polypropylene with polycarbonate window
<b>Filter Type</b>		
Pre-filter	ACFPRE10	ACFPRE10
HEPA	ACFHEPA-188	ACFHEPA-188
<b>Operation</b>		
	Single Speed	Single Speed
<b>Blower</b>		
	Sparkless/Brushless	Sparkless/Brushless
<b>Electrical</b>		
Voltage	120V 60 Hz	220V 50 Hz
Current	0.7A	0.35A
Consumption	84W	84W
<b>Airflow Characteristics</b>		
Volume of Filtered Air	180 CFM	180 CFM



2179 East Lyon Station Road • Creedmoor, NC 27522  
 Phone: (919) 255-9472 • Toll Free: (877) 255-9472 • Fax: (919) 528-3400  
[www.csmedicalllc.com](http://www.csmedicalllc.com)