

TransPorter™

TEE TRANSPORT DEVICE

TransPorter™ Mobile Cart

TransPorter™ Mobile Cart is designed for the simple and safe movement of multiple TPorter™ TEE Ultrasound Probe Transportation and Procedure Cases. TransPorter securely holds four TPorter cases within an individually dedicated compartment. TransPorter mobile cart provides flexibility and a high degree of usability as staff manages the distribution and collection of expensive and delicate TEE ultrasound probes throughout the healthcare facility.

TransPorter has a stainless steel worktop that allows a TPorter case to be placed on top so the technician can position the high-level disinfected TEE ultrasound probe and other critical point-of-care supplies into the dedicated compartments of TPorter. TransPorter is shipped fully assembled with two locking casters and two free-floating casters.



SPECIFICATIONS

CS-200980 TransPorter Mobile Cart	
Width	40"
Depth	25"
Height	40" (work top height = 39")



Creedmoor, NC 27522 USA
 Phone: +1 919-255-9472 • Fax: +1 919-528-3400
www.csmedicalllc.com • sales@csmedicalllc.com