

CleanShield® Ultrasound Probe Storage Cabinet

The CleanShield® Ultrasound Probe Storage Cabinet is designed to effectively and securely store disinfected endocavity transducers and semi-critical devices by suspending each disinfected probe in a vertical position that minimizes stress on the connection cable, strain relief and electrical pack. Specifically designed for vaginal, rectal and general purpose ultrasound probes and for mounting on the wall at point of care.



Cat. # CS-VR06W

SPECIFICATIONS

| | CS-VR06W |
|--------|----------|
| Width | 26" |
| Depth | 10" |
| Height | 60" |
| Weight | 80 lbs. |

SPECIFICATIONS

| | CS-VR06W | CS-VR06W 220 |
|--------------------------------|---|---|
| Construction | | |
| Main Chamber | Polypropylene | Polypropylene |
| Door | Polypropylene with polycarbonate window | Polypropylene with polycarbonate window |
| Filter Type | | |
| Pre-filter | ACFPRE10 | ACFPRE10 |
| HEPA | ACFHEPA-188 | ACFHEPA-188 |
| Operation | | |
| Blower | Sparkless/Brushless | Sparkless/Brushless |
| Electrical | | |
| Voltage | 120V 60 Hz | 220V 50 Hz |
| Current | 0.7A | 0.35A |
| Consumption | 84W | 84W |
| Airflow Characteristics | | |
| Volume of Filtered Air | 180 CFM | 180 CFM |



2179 East Lyon Station Road • Creedmoor, NC 27522
 Phone: (919) 255-9472 • Toll Free: (877) 255-9472 • Fax: (919) 528-3400
www.csmedicalllc.com